


One-Dimensional Barcoding

MED e-care e-MAR now supports one-dimensional barcode functionality. Barcodes can be used to quickly navigate to a resident's current MAR pass.

How to get the barcode?

Navigate to **e-Reports >> ADT >> Group >> Barcode List (Chart Numbers)** and generate the report. This report will display the barcodes for each resident which you can scan and it will navigate to the current MAR pass screen.



Sample Report:

ID 90103

Barcode List

Sector(s): All Sectors
Unit(s): All Units
Program(s): All Programs

Room - Bed	Client Name	Image	Barcode
-	SMITH, Angela		 09EW 13560
-	SMITH, Angelina		 0P1W 36170
-	SMITH, Angelo		 087Y 11852

How it works?

1. Navigate to Settings >> e-MAR >> Options and select "Bar Code Type: Chart Number".
2. Navigate to e-Reports >> ADT >> Group >> Barcode List (Chart Numbers) >> Generate and Print the report.
3. Scan the barcode from the printed report using the provided scanner. (for a list of compatible scanners, have your IT department connect with our support department).
4. Once you scan the code, the system will navigate to the current MAR pass screen where you can administer the medications for the resident in question.
5. If there are no scheduled medications for the resident in the current medication pass, the system displays a warning message.

